

SonoTrax

Ultrasonic Pocket Doppler



SonoTrax
Ultrasonic Pocket Doppler



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.

EDANUSA

4901 Morena Blvd., Suite 505 | San Diego, CA 92117 | Tel: 888.850.4597

EDANUSA
www.edanusa.com

SonoTrax

Ultrasonic Pocket Doppler



SonoTrax II Pro



Compact Battery Charger



Backlit LCD



Built-in Recorder



Interchangeable Waterproof Probe

Specifications

- Display: 45 x 25 mm LCD display
- Heart Rate Range: 50~210 bpm
- Resolution: 1 bpm
- Recording Time: 240 s
- Audio Output Power: 1W
- Auto Power-off: 1 minute with no signal or operation
- Size: 34 (D) x 89 (W) x 141 (H) mm
- Weight: about 300 g (including the battery)



Carry Case



Rolling Stand

- Easy-to-use
- Compact design
- Support 1.5 V AA battery
- Built-in audio recorder
- Large backlit LCD display
- Accurate fetal heart detection with high-fidelity sound
- Support 2 MHz, 3 MHz, 4 MHz, 5 MHz, 8 MHz waterproof probes
- Automatic and manual counting of FHR
- Interchangeable waterproof probe
- Automatic power-off



SonoTrax Vascular



SonoTrax Lite



SonoTrax Basic



SonoTrax Basic A



SonoTrax Pro



SonoTrax II



SonoTrax II Pro